

ZenPure

PureFlo All Nylon Filter Capsules

(SNL Series)

All Nylon Capsule

PureFlo® All Nylon filter capsules are highly retentive, naturally hydrophilic nylon membrane filters that have been specially designed for liquid applications. The nylon membrane in an all-nylon construction provides excellent chemical compatibility and superior flow per unit area as compared to other membrane capsules. No adhesives, binders, or surfactants are used in the manufacturing process. All capsules are rinsed with high-purity water to reduce extractables and downtime.

The naturally hydrophilic nylon membrane does not require pre-wetting with IPA and flushing with DI water, thereby eliminating a potential source of contamination and hazardous waste.

All filter capsules are 100% integrity tested to ensure filter performance each and every time out of the package. PureFlo nylon filter capsules are well-suited for critical applications where superior flow and particle removal efficiency at 0.05, 0.1, 0.2, 0.45, 0.65, 0.8, and 1.2 micron is required.



Applications					
Solvents Water					
Hard Particle	Beverages				
Cell Removal	Pharmaceuticals				
Chemicals Biologic	5				


Specifications


Materials of Construction: Media Supports: None	Media: Naturally Hydrophilic Nylon 6,6				
Cage,	Core, End Caps: Nylon				
	Sealing: Thermally bonded				
Fitting Connections:	See ordering guide for the availability.				
J	(Custom adaptors available upon request)				
Nominal Dimensions:	Lengths: 1.5in. (68mm), 2.5in. (98mm), 5in. (153mm),				
7.5in.	(251mm), and 10in. (281mm)				
Diameter:	2.88 in. (73mm)				
Effective Filtration Area:	1.5 in - 0.74 ft ² (690 cm ²), 2.5 in - 1.4 ft ² (1300 cm ²),				
5.0	in - 2.64 ft ² (2450 cm ²), 7.5 in - 4.0ft ² (3750 cm ²),				
	10 in - 5.3 ft ² (4900 cm ²)				
Operating Conditions:	Maximum Operating Pressure: Liquid: 5.5 bar (80psi) at 72°F/22°C				
Liquid:	3.0 bar (43.5psi) at 203°F/95°C				
Gas:	4.1 bar (60psi) at 72°F/22°C				
	Maximum Forward Differential Pressure: 4.2 bar (60psi) at 72°F/22°C				
	1.5 bar (21.75psi) at 203°F/95°C				
	Maximum Reverse Differential Pressure: 2.1 bar (30 psi) at 72°F/22°C				
	Minimum Burst Pressure: 8.3 bar (120psi) at 72°F/22°C				
Maximum Operating Temperature: : 203°F/95°C					

ZenPure

PureFlo All Nylon Filter Capsules

(SNL Series)






Specifications (cont.)

Autoclavable & Sanitizable:

Capsules can be autoclaved 5 times at 257°F/125°C for 30 minutes or chemically sanitized in situ using common sanitizing agents or hot water at 194°F/90°C for a limited time (dependent on time and temperature). Capsules must not be in situ steam-sterilized.

Regulatory Compliance:

The filters are constructed with polypropylene resins and filtration media in compliance with 21CFR Part 177 of the US Code of Federal Regulations and USP Class VI Biological Test for Plastic.



PureFlo All Nylon Filter Capsule Ordering Guide

PureFlo Capsule Filters	Vent & Drain Fitting	Filter Media	Pore size (Micron)	Length	Input Fitting	Output Fitting	Options
riiteis	L = Female Luer		005= 0.05	H = 1.5"	1C 1C	= 1/8" Compression	Shell Material
SN	Lock	N = Nylon	010 = 0.10	S = 2.5"		1H = 1/8" Hose barbs	Blank = Nylon construction
Media on All Nylon	P = M8 Threaded		020 = 0.20	L = 5"	1N 1N	= 1/8" MNPT	Sterilization
Construction	Plug - Limited		045 = 0.45	E = 7.5"	1NF 1NF		
Omm Diameter	options		065 = 0.65	K = 10"	10 10	= 1/8" Male Quick Coupling with Metal latch	-ETO = Ethylene Oxide Sterilization
			080 = 0.80		1QF 1QF	= 1/8" Female Quick Coupling with Metal latch	Vent & Drain
			120 = 1.20		1QFV 1QF		-N = No vent or Drain fittings
					1QV 1QV	= 1/8" Male Valved Quick Coupling with Metal latch	-NI = No vent or drain inlet fitting
					2H	2H = 1/4" to 1/2" Hose barbs	NO NO CONTRACTOR OUTLAND STATE
					2HS 2HS	= 1/4" to 3/8" Hose barbs	-NO = NO vent or drain Outlet fitting
					2LF	2LF = Large Female Luer Lock fitting	
					2N	2N = 1/4" MNPT	O Diversity Out to Counting
					2NF	2NF = 1/4" FNPT	O-Ring for Quick Coupling
					2P 2P	= 1/4" Push to Connect	Blank = O-ring Silicon (Standard)
					2PE 2PE	= 1/4" Push to Connect Elbow	-OE = O-ring EPDM
					2Q 2Q	= 1/4" Male Quick Coupling for Metal latch	-ON = O-ring Nitrile
					2QF	2QF = 1/4" Female Quick Coupling with Metal latch	-OV = O-ring Viton
					2QFP 2QF	P = 1/4" Female Quick Coupling with Plastic latch	
					2QFV 2QF		
					2QP 2QP		
					2QV 2QV	= 1/4" Male Valved Quick Coupling with Metal latch	
						2V = 1/4" Walther Quick Connect	
					3Н	3H = 3/8" Hose Barb	
					3HE 3HE	= 3/8" Hose Barb Elbow	
					3H-FB		
					3NF 3NF		
						4H = 1/2" Hose Barb	
					4N 4N	= 1/2" MNPT	
					40 40	= 1/2" Male Quick Coupling for Plastic latch	
		I				1/2" Tri clamps	
		I			TC TC	= 1-1/2" Tri clamp	
		I			TS	= Shower Head	
		cample - 0.2um				- Shorter field	

International Filter Products
Revolutionary systems for your industry

150 9001 14001

904 N. Frederic St. Burbank, CA 91505 www.internationalfilterproducts.com sales@internationalfilterproducts.com Phone: (818) 841-2702 Fax: (818) 841-2947

ZenPure and PureFlo are registered trademarks of ZenPure Corporation or an affiliated company. Copyright 2003-2012 ZenPure or an affiliated company. All rights reserved.

