

PureFlo All Nylon Filter Capsules (SNL Series)

All Nylon Capsule

PureFlo® All Nylon filter capsules are highly retentive, naturally hydrophilic nylon membrane filters that have been specially designed for liquid applications. The nylon membrane in an all-nylon construction provides excellent chemical compatibility and superior flow per unit area as compared to other membrane capsules. No adhesives, binders, or surfactants are used in the manufacturing process. All capsules are rinsed with high-purity water to reduce extractables and downtime.

The naturally hydrophilic nylon membrane does not require pre-wetting with IPA and flushing with DI water, thereby eliminating a potential source of contamination and hazardous waste.

All filter capsules are 100% integrity tested to ensure filter performance each and every time out of the package. PureFlo nylon filter capsules are well-suited for critical applications where superior flow and particle removal efficiency at 0.05, 0.1, 0.2, 0.45, 0.65, 0.8, and 1.2 micron is required.





Applications

Solvents	Water
Hard Particle	Beverages
Cell Removal	Pharmaceuticals
Chemicals	Biologics

Specifications

Materials of Construction:	Media: Naturally Hydrophilic Nylon 6,6 Media Supports: None Cage, Core, End Caps: Nylon Sealing: Thermally bonded
Fitting Connections:	See ordering guide for the availability. (Custom adaptors available upon request)
Nominal Dimensions:	Lengths: 1.5in. (68mm), 2.5in. (98mm), 5in. (153mm), 7.5in. (251mm), and 10in. (281mm) Diameter: 2.88 in. (73mm)
Effective Filtration Area:	1.5 in - 0.74 ft ² (690 cm ²), 2.5 in - 1.4 ft ² (1300 cm ²), 5.0 in - 2.64 ft ² (2450 cm ²), 7.5 in - 4.0ft ² (3750 cm ²), 10 in - 5.3 ft ² (4900 cm ²)
Operating Conditions:	Maximum Operating Pressure: Liquid: 5.5 bar (80psi) at 72°F/22°C Liquid: 3.0 bar (43.5psi) at 203°F/95°C Gas: 4.1 bar (60psi) at 72°F/22°C Maximum Forward Differential Pressure: 4.2 bar (60psi) at 72°F/22°C 1.5 bar (21.75psi) at 203°F/95°C Maximum Reverse Differential Pressure: 2.1 bar (30 psi) at 72°F/22°C Minimum Burst Pressure: 8.3 bar (120psi) at 72°F/22°C Maximum Operating Temperature: : 203°F/95°C

PureFlo All Nylon Filter Capsules (SNL Series)




Specifications (cont.)

Autoclavable & Sanitizable:

Capsules can be autoclaved 5 times at 257°F/125°C for 30 minutes or chemically sanitized in situ using common sanitizing agents or hot water at 194°F/90°C for a limited time (dependent on time and temperature). Capsules must not be in situ steam-sterilized.

Regulatory Compliance:

The filters are constructed with polypropylene resins and filtration media in compliance with 21CFR Part 177 of the US Code of Federal Regulations and USP Class VI Biological Test for Plastic.



PureFlo All Nylon Filter Capsule Ordering Guide

PureFlo Capsule Filters	Vent & Drain Fitting	Filter Media	Pore size (Micron)	Length	Input Fitting	Output Fitting	Options
SN Media on All Nylon Construction 70mm Diameter	L = Female Luer Lock P = M8 Threaded Plug - <u>Limited options</u>	N = Nylon	005 = 0.05 010 = 0.10 020 = 0.20 045 = 0.45 065 = 0.65 080 = 0.80 120 = 1.20	H = 1.5" S = 2.5" L = 5" E = 7.5" K = 10"	1C 1H 1N 1NF 1Q 1QF 1QFV 1QV 2H 2HS 2LF 2N 2NF 2P 2PE 2Q 2QF 2QFP 2QFV 2QP 2QV 2V 3H 3HE 3NF 4H 4N 4Q MT TC	1C = 1/8" Compression 1H = 1/8" Hose barbs 1N = 1/8" MNPT 1NF = 1/8" FNPT 1Q = 1/8" Male Quick Coupling with Metal latch 1QF = 1/8" Female Quick Coupling with Metal latch 1QFV = 1/8" Female Valved Quick Coupling with Metal latch 1QV = 1/8" Male Valved Quick Coupling with Metal latch 2H = 1/4" to 1/2" Hose barbs 2HS = 1/4" to 3/8" Hose barbs 2LF = Large Female Luer Lock fitting 2N = 1/4" MNPT 2NF = 1/4" FNPT 2P = 1/4" Push to Connect 2PE = 1/4" Push to Connect Elbow 2Q = 1/4" Male Quick Coupling for Metal latch 2QF = 1/4" Female Quick Coupling with Metal latch 2QFP = 1/4" Female Quick Coupling with Plastic latch 2QFV = 1/4" Female Valved Quick Coupling with Metal latch 2QP = 1/4" Male Quick Coupling for Plastic latch 2QV = 1/4" Male Valved Quick Coupling with Metal latch 2V = 1/4" Walther Quick Connect 3H = 3/8" Hose Barb 3HE = 3/8" Hose Barb Elbow 3H-FB = 3/8" Hose Barb with filling bell or -FC with Cap 3NF = 3/8" FNPT 4H = 1/2" Hose Barb 4N = 1/2" MNPT 4Q = 1/2" Male Quick Coupling for Plastic latch MT = 1/2" Tri clamps TC = 1-1/2" Tri clamp TS = Shower Head	Shell Material Blank = Nylon construction Sterilization -ETO = Ethylene Oxide Sterilization Vent & Drain -N = No vent or Drain fittings -NI = No vent or drain inlet fitting -NO = NO vent or drain Outlet fitting O-Ring for Quick Coupling Blank = O-ring Silicon (Standard) -OE = O-ring EPDM -ON = O-ring Nitrile -OV = O-ring Viton
Example - 0.2um Nylon with Nylon Construction, Luer Lock V/D, 1.5" internal filter, 1/4" MNPT I/O is SNLN020H2N2N							
Quick Couplings are molded, Female Quick Coupling with Plastic latch will be threaded on, compatible with CPC (Colder), LinkTech and others.							

Your Local Distributor:

ZenPure

North & South Americas:
ZenPure Americas, Inc
www.zenpure.com
Info-us@zenpure.com
703-335-9910

All Other Regions:
ZenPure Corporation
www.zenpure.com
info@zenpure.com
+86 571 2288 6800

ZenPure and PureFlo are registered trademarks of ZenPure Corporation or an affiliated company. Copyright 2003-2012 ZenPure or an affiliated company. All rights reserved.

