

ZenPure

PureFlo® Glass Microfiber Filter Capsules

PureFlo® Borosilicate Microfiber filter cartridges are highly retentive filters that have a pleated design to maximize surface area. The borosilicate microfiber media in an all-polypropylene construction provides excellent chemical compatibility and superior flow per unit area. No adhesives or surfactants are used in the manufacturing process. The borosilicate microfiber is "non-fiber releasing."

The borosilicate microfiber is ideally-suited for aggressive microbial retention. In addition, the media provides superior flow and pressure drop characteristics per unit area. PureFlo Microfiber filter cartridges are well-suited for critical applications where superior flow and particle retention for deformable and non-deformable particle removal is required.

Applications						
Prefiltration	Biologics					
Serums	Bio-burden reduction					
Pharmaceuticals	Retention of deformable particles					
Wine Clarification	Inks					
Food & Beverage	Cosmetics					

Regulatory Compliance

Manufactured from materials that conform to the requirements of USP Class VI Biological Test for Plastics.



Materials of Construction

Media: Non-Woven Borosilicate Microfiber Pleate Supports: Polypropylene Cage, Core, End Caps: Polypropylene Fitting Connections: NPT, Couplings Compression, Hose Barb, Sanitary Flange

Dimensions (TC fittings)

Length: 1.5 in (128 mm) Length: 2.5 in (155 mm) Length: 5.0 in (210 mm) Length: 7.5 in (270 mm)

Length: 10 in (340mm)
Diameter: 2.75 in (70 mm) **Effective Filtration Area**

1.5 in - 0.48 ft² (450 cm²) 2.5 in - 1.0 ft² (950 cm²) 5.0 in - 1.8 ft² (1700 cm²) 7.5 in - 4.3 ft² (3980 cm²)

10 in - $3.66 \text{ ft}^2 (3400 \text{ cm}^2)$

Operating Conditions

Maximum Operating Pressure:

Liquid: 5.5 bar (80psi) at 77°F/25°C Gas: 4.1 bar (60psi) at 77°F/25°C

Minimum Burst Pressure:

8.3 bar (120psi) at 77°F/25°C

Maximum Forward Differential Pressure:

4.21 bar (60psi) at 68°F/20°C

Maximum Reverse Differential Pressure:

2.1 bar (30psi) at 68°F/20°C

Maximum Operating Temperature:

PP & Gamma PP: 176°F/80°C

HDPE: 140°F/60°C

Features	Benefits
■ Borosilicate Microfiber Media (Absolute Rated)	 High dirt-loading capacity Natural positive charge Superior process throughputs
■ Wide Chemical & Thermal Compatibility	 Provides excellent compatibility No media migration
■ Quality Construction	 No adhesives or surfactants Thermally welded construction

ZenPure

PureFlo® Glass Microfiber Filter Capsules



Inlet/Outlet Fittings	Vent and Drain fittings						
1H = 1/8" Hose barbs	Luer Lock						
1NF = 1/8" FNPT	Luer Lock						
2C = 1/4" Compression	1/4" Comp						
2H = 1/4" to 1/2" Hose barbs	Luer Lock						
2HS = 1/4" to 3/8" Hose barbs	Luer Lock						
2L = Large Female Luer Lock fitting	Luer Lock						
2N = 1/4" MNPT	Luer Lock						
2NF = 1/4" FNPT	Luer Lock						
2Q = 1/4" Male Quick Coupling	Luer Lock						
2QF = 1/4" Female Quick Coupling	Luer Lock						
2QFP = 1/4" Female Quick Coupling with Plastic latch	Luer Lock						
2QP = 1/4" Male Quick Coupling for Plastic latch	Luer Lock						
3C = 3/8" Compression	1/4" Comp						
3H = 3/8" Hose Barb (future)	Luer Lock						
3HE = 3/8" Hose Barb Elbow	Luer Lock						
3NF = 3/8" FNPT	Luer Lock						
4C = 1/2" Commpression	1/4" Comp						
4H = 1/2" Hose Barb (future)	Luer Lock						
4N = 1/2" MNPT	Luer Lock						
4QP = 1/2" Male Quick Coupling for Plastic latch	Luer Lock						
MT = 0.5" Tri clamps	Luer Lock						
TC = 1-1/2" Tri clamp	Luer Lock						
Option: Plugged and Cut off Vent and/or Drains							

PureFlo Glass Fiber Filter Cartridge Ordering Guide

PureFlo Capsule Filters	Vent & Drain Fitting	Filter Media	Pore size (Micron)	Length	Input Fitting	Output Fitting	Options	
sĸ	L = Female Luer Lock	- 0		_	1H = 1/8" Hose barbs 1NF = 1/8" FNPT	1H = 1/8" Hose barbs 1NF = 1/8" FNPT	-GP = Gamma stable polypropylene shell	O-Ring for Quick Coupling
70mm Diameter Capsule	Dec. Person of the control of the co	G = Glass Fiber	010 = 1.0 030 = 3.0	L = 5" E= 7.5" K = 10"	1Q = 1/8" Male Quick Coupling for Metal latch 1QF = 1/8" Female Quick Coupling w/Metal latch 2QF = 1/4" Compression 2H = 1/4" to 1/2" Hose barbs 2HS = 1/4" to 1/2" Hose barbs 2HS = 1/4" to 3/8" Hose barbs 2HF = Large Female Luer Lock fitting 2N = 1/4" MNPT 2NF = 1/4" FNPT 2QF = 1/4" FNPT 2QF = 1/4" Female Quick Coupling my/Metal latch 2QF = 1/4" Female Quick Coupling w/Plastic latch 2QF = 1/4" Hale Quick Coupling my/Plastic latch 2QF = 1/4" Male Quick Coupling for Plastic latch 3QF = 1/8" Female Quick Coupling for Plastic latch 3G = 3/8" Compression 3H = 3/8" Hose Barb 3HE = 3/8" Hose Barb Elbow 3NF = 3/8" FNPT 4C = 1/2" Compression 4H = 1/2" Hose Barb 4N = 1/2" MNPT	1Q = 1/8" Male Quick Coupling for Metal latch 1QF = 1/8" Female Quick Coupling w/Metal latch 2QF = 1/4" Compression 2H = 1/4" to 1/2" Hose barbs 2HS = 1/4" to 3/8" Hose barbs 2HS = 1/4" to 3/8" Hose barbs 2LF = Large Female Luer Lock fitting 2N = 1/4" NNPT 2NF = 1/4" FNPT 2QF = 1/4" FPME 2QF = 1/4" Female Quick Coupling w/Metal latch 2QF = 1/4" Female Quick Coupling w/Plastic latch 2QF = 1/4" Female Quick Coupling w/Plastic latch 2QF = 1/4" Female Quick Coupling for Plastic latch 3QF = 3/8" Hose Barb 3HE = 3/8" Hose Barb Elbow 3NF = 3/8" FNPT 4C = 1/2" Compression 4H = 1/2" Hose Barb 4N = 1/2" MNPT	-NY = Nylon Construction -FB = Filling bell -FC = Filling bell w/Cap -ETO = Ethylene oxide sterilization -NI = No vent or drain inlet fitting -NO = NO vent or drain outlet fitting -E = Polyethylene shell	-OE = O-ring EPDM -ON = O-ring Nitrile -OV = O-ring Viton Blank = O-ring Silicon (Standard) PreFilters -G(pore Size) = Glass Fiber PreFilter -P(pore Size) = PolyPro Media PreFilter -S(pore Size) = PES PreFilter
		Example - SK	Series Capsule,		4QP = 1/2" Male Quick Coupling for Plastic latch MT = 1/2" Tri clamps TC = 1-1/2" Tri clamp //D, 2.5" Internal filter. 0.3 micron Glass Fiber Capsul	4QP = 1/2" Male Quick Coupling for Plastic latch MT = 1/2" Tri clamps TC = 1-1/2" Tri clamp 2, Standard Grade, with 1.5" Tri-Clamp fittings would be	SKLG030STCTC	

International Filter Products
Revolutionary systems for your industry

904 N. Frederic St. Burbank, CA 91505 www.internationalfilterproducts.com sales@internationalfilterproducts.com Phone: (818) 841-2702

Fax: (818) 841-2947

ZenPure and PureFlo are registered trademarks of ZenPure Corporation or an affiliated company. Copyright 2003-2012 ZenPure or an affiliated company. All rights reserved.



